nanoleaf®
Smarter by Design


6-in-1 Light Therapy Wand

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Quick Start Guide
FR Guide de démarrage rapide **DE** Schnellstartanleitung **IT** Guida rapida all'uso **ES** Guía de inicio rápido **PT** Guia de início rápido **BG** Ръководство за бърз старт **CS** Průvodce rychlým startem **DA** Hurtigstartvejledning **EL** Οδηγός γρήγορης εκκίνησης **ET** Kiirjuhend **FI** Pika-aloitusopas **HR** Brzi vodič za početak **HU** Gyors üzembe helyezési útmutató **LT** Greito paleidimo vadovas **LV** Ātrās sākšanas ceļvedis **NL** Snelstartgids **NO** Hurtigstartveiledning **PL** Przewodnik szybkiego startu **RO** Ghid de pornire rapidă **SK** Stručný sprievodca spustením **SV** Snabbstartsguide **SC** 快速入门指南 **TC** 快速入門指南 **KR** 빠른 시작 가이드 **JP** クイックスタートガイド

What's in the Box?

Contenu de la boite

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
6 modes de Traitement

Anti-Aging

Red Light + Near Infrared (NIR) + Heating + Vibration

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Soothe
Heating + Vibration
3 Intensity Levels

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*It is normal for the wand's metal head to be slightly warm after use of the Cooling Mode.
After using the Repair or cooling modes, you need to wait 4 minutes before using it again.

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
Important Things to Note

Do not use the device while it is charging.
The 3 LEDs above the Intensity Button also help to indicate the battery level.

Fully Charged


All 3 LEDs stay on when plugged in and fully charged. Takes 4hrs for a full charge.

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Charging


The 3 LEDs cycle from bottom to top. It cannot be used while charging.

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Low Battery

The 3 LEDs cycle from top to bottom, and beep every 5s.

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
À savoir

WARNING

This is not a medical device and cannot be used for professional treatment or as a substitute for professional medical advice. The actual effect of use varies from person to person.

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


Mode d'emploi

Wash your face before use.

Optional: Apply your favorite serums and moisturizers.

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Suggested Wand Movements

Target specific areas of your face, or follow any of the suggested wand movements below.

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EN WARNINGS Do not leave the device unattended when turned on to avoid the risk of fire or burns; This device is not designed for use by persons with limited physical, sensory, or mental abilities, or those lacking experience and knowledge, unless accompanied or supervised by someone to avoid fire hazards or burns. This device is not intended for use by children. Children should be supervised to ensure that they do not play with the equipment to avoid the risk of fire and burns. Do NOT wash electrical parts with water or other liquids. Water and electricity are a dangerous combination, and in order to avoid the danger of the exit point, Do not use the device in a humid environment (e.g., in a bathroom or shower or pool attachment) Do not allow water to enter the machine. If the main cord is damaged, you shall have it replaced by Nanogrid Limited, a service centre authorized by Nanogrid Limited or similarly qualified persons in order to avoid risk of electrocution. Do not look directly at the light-emitting diodes when using this device. Do not use the device if it appears to be operating incorrectly, and contact customer support. Do not use the device in combination with muscle heating cremes or other medications which may result in damage to skin tissue. Do not expose any parts of the device to chemical solvents or harsh cleaning fluids. The device contains silicon and plastics that come into contact with the skin during use. Do not use the device if you have a known allergic reaction to contact with silicone or plastics. Avoid fluids on the light panel and controller. Discontinue use if there are water in the device. Opening or altering the device may damage the operating circuits and will result in voiding the warranty. Please keep the unit out of reach of children, since the strap can cause strangulation. Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. Use of controls or adjustments or performance of procedures other than those specified herein may result in HAZARDOUS radiation exposure. Do not store or use this product in high temperature. Do not subject this product to physical damage from external impact or a high fall. The controller of this device contains lithium batteries, and personnel without adequate training should not replace the batteries without authorization to avoid related risks. This device contains a Lithium-ion Battery. Before disposal of this device, please review the local laws and requirements related to Lithium-ion Battery disposal at your area before disposing the device. The preferred method of disposal is recycling. Do not subject the device to physical shocks, impacts, or drop Inspect before each use:Check for dents, cracks, or deformation. Do not use if structural integrity is compromised.

WARNINGS: SAFE USE OF BATTERIES This product uses lithium batteries. Do not place or store this product in a humid or high temperature environment for a long time. Do not dispose of battery into fire or dismantle it under any condition. Do not attempt to disassemble or modify the battery. Doing so can lead to leakage, fire, or explosion, and will void the warranty. Proper disposal and recycling methods should be followed for used batteries. Always refer to the manufacturer's guidelines for battery care and maintenance.
Do not leave your device's battery unused for extended periods. For optimal battery health, charge the battery to 50% at least once every three months to maintain your battery performance.

CE INFORMATION NanoGrid Limited (Nanoleaf®) hereby declares that this device is in compliance with the essential requirements and other relevant provisions of Directive 2014/30/EU. The declaration of conformity is available at https://nanoleaf.me/eudocs. **CLEANING AND MAINTENANCE** Do not clean the device with liquid or sprays when it is connected to power. Clean with a soft dry cloth, and do not use harsh cleaning agents.

WARRANTY INFORMATION This product has a one year warranty from the date of purchase. Damages caused by time and use is not within the scope of the service guarantee. The warranty scope of this product does not include failures cased by any factors other than the defects of the product itself.

Disposal of your old product: Your product is designed and manufactured with high-quality materials and components, which can be recycled and reused. When this crossed-out wheelie bin symbol is attached to a product, it means that the product is covered by European Directive 2012/19/EU. Please find out about your local collection system for electrical and electronic products. Please act according to your local rules and do not dispose of your old products with your normal household waste. The correct disposal of your old product will help to prevent potential negative consequences for the environment and human health.

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Nanoleaf Europe SAS