nanoleaf

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Controls

Power Button

Hold for 2s to Turn On /Off
By default, your device is set to turn off after 10 minutes.
Press the power button once more when on to pause the timer.

Mode Selection

Press once to switch between Red, Green, Blue, Cyan, Yellow, Purple, & White light therapy.

Each light mode is combined with Near Infrared Light (NIR) to enhance the treatment effects. The screen will display the corresponding light mode.
Double-press to cycle between 5 different levels of intensity.
The screen will display the corresponding light intensity.

Timer

Press once to add 1 minute to your timer.
Long press to add minutes faster, up to a max of 30 minutes.
Your device will automatically turn off when your timer is up.

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
Light Therapy Modes

Choose the ideal treatment mode based on your skin's needs.

Red Light (Mode R)
640nm + 850nm NIR
+ 940nm NIR + 1072nm NIR
Visibly smooths wrinkles for a firmer skin appearance.

Green Light (Mode G)
525nm + 850nm NIR + 940nm NIR
Supports a more even-looking skin tone.

Blue Light (Mode B)
460nm + 850nm NIR + 940nm NIR
Targets visible blemishes, leaving skin clearer and refined.

Cyan Light (Mode C)
460nm + 525nm + 850nm NIR + 940nm NIR
Soothes stressed-looking skin for a calm, balanced complexion.

Yellow Light (Mode Y)
525nm + 640nm + 850nm NIR
+ 940nm NIR + 1072nm NIR
Softens visible redness for a more even, refined tone.

Purple Light (Mode P)
460nm + 640nm + 850nm NIR
+ 940nm NIR + 1072nm NIR
Refines rough, uneven texture for a smoother, softer skin surface.

White Light (Mode L)
460nm + 525nm + 640nm + 850nm NIR
+ 940nm NIR + 1072nm NIR
Enhances firmness and softens visible signs of aging.

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
1072nm Deep NIR Treatment Zones


1072nm Deep NIR is strategically positioned to deliver focused care to key areas of concern.

Red Light
(Mode R)

Yellow Light
(Mode Y)

Purple Light
(Mode P)

White Light
(Mode L)

The dedicated 1072nm Deep NIR lights activate only in Red, Yellow, Purple, and White Modes to deliver focused care to key areas of the neck and décolletage. They remain OFF in Blue, Cyan, and Green Modes, which are designed to gently soothe the skin. **This is intentional and does not indicate a product defect.**

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


How to Attach the Straps
Comment fixer les sangles

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
Things You Need to Know
Informations importantes

Do not use the device while it is charging.

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WARNING
AVERTISSEMENT

This is not a medical device and cannot be used for professional treatment or as a substitute for professional medical advice. The actual effect of use varies from person to person.

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

Plug the Neck Mask into the controller and press the power button to turn it on. This product is rechargeable.

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

https://go.nanoleaf.me/myneckmask

FR **Besoin d'aide?** https://go.nanoleaf.me/myneckmask

Nanoleaf Canada
#200–106 Front Street East
Toronto, ON
M5A 1E1, Canada

Nanoleaf Europe SAS
11 Rue de Lourmel
75015 Paris
France
Hello@nanoleaf.me

Nanoleaf Europe SAS
Ruby House, 8 Ruby Place
Aberdeen AB10 1ZP
Great Britain

Nanoleaf Hong Kong
#1405 - 135 Bonham
Strand Trade Centre,
Hong Kong

Model: E50
Rated Voltage: DC 5V
Rated Power: 5W

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