nanoleaf

Lymphatic Body Sculpting Massager
Appareil de massage lymphatique sculptant

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
Controls

Power On / Off
Press and hold the Power Button for 2 seconds.

Select Intensity
Short press the corresponding function button to increase its level.
Quickly press twice to jump 5 levels.
Red light is controlled through the Heat button.

Low Battery
When the battery level is below 20%, the white light will flash as a reminder to charge the device.

Fully Charged
The white LED remains steadily lit when the device is fully charged.


Light Care | Levels 1–10 | 38°C
For first-time use or more gentle.

Routine Care | Levels 11–20 | 40°C
For regular body care.

Deep Care | Levels 21–30 | 42°C
For a stronger body care session.

Choose the level that feels right for your body and comfort. The device automatically turns off after 30 minutes.

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Microcurrent Sculpting

Gentle microcurrent pulses work to firm and sculpt the body's contour.

Vibration Massage

Rhythmic massage encourages lymphatic flow for a smoother, refreshed look.

Heat + Advanced Light Therapy

Comforting warmth pairs with 640nm red light and 850nm near-infrared light to help revitalize skin and promote a smoother, firmer appearance.

Sculpting Mode may feel stronger for some users.
Always begin at Level 1 and increase gradually to a comfortable level. More intensity is not always more effective.

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Even Pressure Distribution

Delivers balanced pressure across the body's curves.

Targeted Stimulation Nodes

Textured nodes focus on specific areas for more precise massage.

Contour-Adaptive Design

Maintains consistent contact across curved body areas.

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

2

Select Mode


Use Sculpting Mode on the abdomen, thighs, and upper arms. Use Vibration Mode for shoulders or neck.

3

Glide Gently


Move slowly with light pressure. Do not stay on one spot for too long in Sculpting Mode. (Refer to “Where to Use” for direction.)

4

Clean After Use

Wipe the device head with a soft, dry cloth after use.

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







Thigh
Glide from the upper thigh downward.

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Where to Use
Où utiliser


Abdomen

Glide downward, or upward for sagging concerns, and use circular motions.

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
Legs/Feet

Legs: Glide upward from the ankle toward the knee.
Feet: Massage from the heel toward the toes.

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
Important Things to Note

Start low and increase gradually

Always choose the care level that feels right for your body and comfort.

Temporary redness or mild tenderness may occur

If bruising appears, lower the level, reduce pressure, and allow the area to rest before using again.

Use Sculpting Mode on body-contouring areas only

DO NOT use Sculpting Mode on shoulders, neck, joints, or bony areas. Use Vibration Mode for shoulders or neck.

Do not press too hard

More pressure does not mean better results. Use light pressure and keep the device moving.

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

Important Things to Note
Points importants à noter


Charge
Recharge

3 hrs for full charge.
This product is rechargeable.

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Portable
Portable

If the device becomes excessively hot, emits smoke, or produces an unusual odor during charging, please immediately disconnect the charger and contact customer service.

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
Important Things to Note
Points importants à noter

WARNING
AVERTISSEMENT

This is not a medical device and cannot be used for professional treatment or as a substitute for professional medical advice. The actual effect of use varies from person to person.

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Avoid looking directly at the LED lights.

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

Do not drop or throw the product on hard surfaces.

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Nanoleaf Canada
#200–106 Front Street East
Toronto, ON
M5A 1E1, Canada

Nanoleaf Europe SAS
11 Rue de Lourmel
75015 Paris
France
Hello@nanoleaf.me

Nanoleaf Europe SAS
Ruby House, 8 Ruby Place
Aberdeen AB10 1ZP
Great Britain

Nanoleaf Hong Kong
#1405 – 135 Bonham
Strand Trade Centre,
Hong Kong

Model: XY-406
Power Rating: 6W
Input: DC 5V
Battery Capacity: 1800mAh

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